



EEG & fNIRS Based Alzheimer's Disease Detection

Primarily references the literature: *Event-specific EEG-FNIRS feature fusion FOR Alzheimer's disease classification*^[1]

Gangfeng Hu, Sangdae Nam, Niko Hams,
Uli Prantz, Eric Ji, Matthew Zhou

Introduction: Alzheimer's Disease

- An irreversible neurodegenerative disease
- Starts slowly and progressively worsens
- 60–70% of cases of all dementia cases.




Exact cause is not fully understood:

Some **key factors** believed to be involved:


- Genetics
- Abnormal protein deposition: amyloid- β and tau
- Vascular factors
- ...

Introduction: Alzheimer's Disease


- **AD** (Alzheimer's Disease): Severe memory and cognitive decline
- **MCI** (Mild Cognitive Impairment): Noticeable but mild cognitive changes
- **NC** (Normal Cognition): Healthy brain function, no impairments
- Focus of Study: Classify participants into these three categories



Normal




MCI



AD

Introduction: EEG & fNIRS



Introduction: EEG & fNIRS


Why Combine EEG and fNIRS?

- Exploit complementary strengths

EEG	fNIRS
Real-time	Slight delay in response
High temporal	Less motion-sensitive
Less spatial resolution	High spatial resolution



EEG & fNIRS



Methods: the Four Tasks


EEG and fNIRS signals are detected during the process of the four tasks below:

1. Resting Task:

- Focus on a white cross on a screen
- Duration: 60 seconds
- Baseline brain activity measurement


3. 1-back Task:

- A number in {1,2,3} shows on a screen every time
- Press a button when the current number matches the previous one
- Evaluates working memory and cognitive flexibility



2. Oddball Task:


- Identify yellow (target) and blue (non-target) circles
- Respond to yellow by pressing a button
- Tests attention and response inhibition



4. Verbal Fluency Task:


- **Phonemic Fluency:** Generate words starting with a given letter
- **Semantic Fluency:** Generate words within a specific category (e.g., animals)

Methods: Pipeline



Methods: Feature Extraction


(a) EEG signal feature extraction



PSD (Power Spectral Density, W/Hz or dB/Hz): for a certain segmented signal $x(t)$, let $X(f)$ be its Fourier transform: then the PSD of the signal is

$$P(f) = |X(f)|^2$$

(b) fNIRS signal feature extraction



Experiment Results

Performance Metrics:

Metric	ExtraTree	MLP
Accuracy	0.79	0.83
Macro Avg F1	0.80	0.83
Weighted Avg F1	0.79	0.83

Key Observations

- MLP performs slightly better in overall accuracy and F1 scores.
- ExtraTree excels in precision for Class 0 but has lower recall for Class 2.
- MLP balances recall better across classes, leading to higher overall performance.